




サルファ剤
Sulfonamides 磺胺药物


O950306E




Sulfaguanidine
磺胺胍




Sulfapyridine
磺胺吡啶




Sulfanilamide
磺胺




Sulfamerazine
磺胺甲噁啶




Sulfisomidine
磺胺索嘧啶




Sulfamethizole
磺胺甲唑




Sulfathiazole
磺胺塞唑



Sulfamethoxy pyridazine
磺胺甲氧哒嗪



Sulfadiazine
磺胺嘧啶



Sulfamethazine
磺胺二甲嘧啶

Column : J'sphere ODS-M80 (4 μm, 8nm)
250 × 4.6 mm I.D.
Eluent : THF/50mM NH₄H₂PO₄ (10/90)
Flow rate : 0.5 mL/min
Temperature : 37 °C
Detection : UV at 275 nm, 0.16 AUFS
Injection : 3 μL (50 μg/mL)